

One Cell Multiple Readouts with FocuSCOPE®


High-throughput single-cell sequencing technologies have significantly enhanced our understanding of complex biological systems. However, most current high-throughput single-cell sequencing methods detect only gene expression levels. Mutation detection at single-cell resolution remains a bottleneck. To address this, Singleron has developed FocuSCOPE, a high-throughput multi-omics sequencing solution capable of **detecting mutations, gene fusions, viral sequences, and the whole transcriptome** from the same single cells, better elucidating the connection between genotype and gene expression.

Highlights

- * High sensitivity
- * High accuracy
- * Not limited to poly-T capture

Principle

The core technology of FocuSCOPE® is centered around specially designed **Barcoding Beads** containing customizable probes and poly T probes that capture poly adenylated mRNA and specific transcripts. During the workflow, cDNA is split, and two different library is prepared which can be combined for sequencing.



4 ready to use panel with a possibility to fully customize!


Lung cancer*	Clonal hematopoiesis*	Blood cancer*	Epstein-Barr Virus*	Custom
* EGFR	* DNMT3A	* WT1	* EBNA1	
* KRAS	* TET2	* KRAS	* EBNA2	
* PIK3CA	* ASXL1	* IDH1/IDH2	* EBER1	
* BRAF	* JAK2	* TP53	* EBER2	
* TP53	* TP53	* BCR_ABL1	* ZEBRA	
		* PML_RARA		

*Specific mutation sites are targeted in the genes showed above.

FocuSCOPE exhibits high degree of specificity


NB4, CCRF and K562 were mixed in equal proportions. Libraries were prepared by FocuSCOPE® Single Cell Multiomics Blood Cancer mRNA x Mutation Analysis Kit.

NB4 cell line contains KRAS (A18D), **CCRF** cell line contains KRAS (G12D) mutations.



High capture rate compared to poly-T capture

NB4 cell line that contains PML-RARA fusion gene was subcutaneously injected into immunodeficient mice to generate CDX (cell-line-derived xenograft) tumor model. Libraries are generated by using FocuSCOPE BC Barcoding Beads or poly-T beads.



Bead type	Nr. NB4 cells	Nr. of PML-RARA	Detection rate
Poly-T Beads	11504	724	6%
BC Barcoding Bead	11543	6876	59%

Contact us to get more information!

Singleron Biotechnologies GmbH

☎ Tel.: +49 221 16824777 ✉ Email: info@singleronbio.com 📍 Address: Gottfried-Hagen-Str. 60, 51105 Cologne, Germany